

# Flat Digital Scale



Accurate results of the DS600 are easy-to-read thanks to a clear display and 2" tall numbers.



## Sleek design with easy-to-read display.


### KEY FEATURES

The Auto-On feature of the DS600 switches on the scale when the patient taps on the heavy duty glass platform. The weight is shown in a 2" high LCD display. The tempered glass weighing surface is 12" x 12" and the height of the scale is just 3/4" high. The DS600 can display weight in lb or kg. The DS600 is EMR compliant.

### iWELLNESS APP

Track your health and stay fit with the new iWellness app. Easily connect the app to your DS600 using Bluetooth. The app integrates with the iOS Health app as a data source for weight. All your data will automatically be recorded into the app.

### SPECIFICATIONS

<b>Warranty:</b>	Two Years
<b>Maximum Capacity:</b>	400 x 0.2 lb (180 x 0.1 kg)
<b>Functions:</b>	Auto On / Off, BMI
<b>Communication:</b>	Bluetooth 4.0 – iOS and Android compatible
	
<b>Display:</b>	2" High LCD
<b>Construction:</b>	6 mm Tempered Glass
<b>Dimensions:</b>	12.75 L x 12.75" W x 1" H
<b>Displayed Units:</b>	lb, kg, St
<b>Power:</b>	4 AA Batteries (Included)

Download iWellness for iPhone or Android

