

Flat Digital Scale

KEY FEATURES

The *Auto-On* feature of our DS500 switches the scale “on” when the patient taps on the heavy duty glass platform. A patient’s accurate weight measurement is easy-to-read because of the unit’s 2” high back-lit LCD display. The durable tempered glass weighing surface is 12” x 12” and the scale’s sleek contemporary design is just 3/4” high. The DS500 can display weight in lbs or kg.

- Auto-On function for easy start up
- Spacious weighing platform
- Convenient 0.75” scale height
- Durable tempered glass construction
- Easy-to-read LCD display with 2” numerals
- Quickly switch between lb or kg

SPECIFICATIONS

Warranty:	Two Years
Power Input:	3V Lithium Battery (Included)
Battery Life:	3,000 Weighments
Maximum Capacity:	400 x 0.2 lb (180kg x 0.1kg)
Accuracy:	+/- 0.2 lb (0.1 kg)
Construction:	Heavy Duty Tempered Glass
Dimensions:	11.75" L x 12.25" W x 0.75" H
Displayed Units:	lb, kg

