

# Handrail Scale

## KEY FEATURES

Accuracy, patient safety, and portability have all been engineered into the Doran DS7100 Handrail Scale. Because the scale only weighs 30 pounds, moving it from room to room takes very little effort. Simply tilt the scale back on the two built-in, non-marking wheels and roll to your destination.

As soon as the patient steps onto our well designed, low-profile DS7100 platform, an accurate weight up to 1,000 pounds will be held on the display. Doran's exclusive *Weight Lock* feature ensures an accurate weight is displayed regardless of a patient's movements. Our *Reweigh* functionality allows a second weight to be acquired while the patient remains on the platform.

After the patient has stepped off the platform, *Last Weight* allows the previous weight and BMI value to be recalled for later recording.



### DS7100-HR

The Model 7100-HR offers the same features of the 7100 Handrail Scale with the addition of an integrated height rod.

## SPECIFICATIONS

<b>Warranty:</b>	Three Years
<b>Power Input:</b>	Six AA Batteries (Included)
<b>Battery Life:</b>	11,000 Weighments
<b>Maximum Capacity:</b>	1000 x 0.1 lb (454 x 0.05 kg)
<b>Accuracy:</b>	+/- 0.1 lb (0.05 kg)
<b>Construction:</b>	Aluminum and Painted Steel with Plastic Trim
<b>Dimensions:</b>	Scale Dimensions 25" W x 29" D x 46" H  Platform Dimensions 22.5" W x 20.5" D x 2" H
<b>Display:</b>	1" High LCD
<b>Displayed Units:</b>	lb, kg, Units lock feature
<b>Controls:</b>	Zero, LB, KG, Last Weight, Reweigh, Scroll Up/Down and BMI
<b>BMI Calculator:</b>	Input data in FT/IN or CM
<b>Options:</b>	RS-232 Communications Port USB Communications Port AC Adapter International Power Supply

