

# Patient Weight Data Collection

## Excellerator Software



Manually recording patient weight information and later retyping this data into a spreadsheet is now a thing of the past. Automate your data collection and remove the risk of human error with the easy to use Excellerator Data Collection Software. Collecting patients' weight and BMI data from the scale is quick and easy with the optional USB or RS-232 interface. Connect the scale to a PC with Microsoft Excel installed and spreadsheets will be populated with data in real-time. Once configured, you can begin collecting data into Microsoft Excel in just two clicks of the mouse, it's really that easy.



Connect a Doran Medical Scale to a PC with the available USB Communications Option.

### MINIMUM SYSTEM REQUIREMENTS

- 133 MHz or faster processor
- 64 Mb of RAM
- CD-ROM Drive
- Minimum 20Mb of available hard drive space
- Windows 2000 or Windows XP Operating System
- Microsoft Excel Version 97 or higher
- One available RS232 serial or USB communications port



**Doran Scales, Inc.** 1-800-264-4106  
1315 Paramount Parkway Fax: 630-879-0073  
Batavia, IL 60510 [www.doranmedical.com](http://www.doranmedical.com)